

Urine Bag

No Outlet with Hanger

DESCRIPTION

Our basic Urine Bags are used to collect and contain urine via a catheter or leg bag and are made from durable material. They also feature clear volume graduation to aid simple volume measurement. They are available with hanger.

FEATURES & BENEFITS

- Durable material
- Clear volume graduation
- With non-return valve
- With hanger
- No outlet


PRODUCT DATA

- Tube length: 120cm
- Tube ID: 6.5mm
- Tube OD: 9mm
- 2000mL Bag
- Material: PVC, PP, PE
- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Shelf life: 5 years
- ARTG Number (no outlet): 271114
- Sterile EO

ISO STANDARDS

- ISO 14971:2012 Medical Device - Applications of Risk Management.
- ISO 10993-1:2009 Biological compatibility test report
- ISO 8669:1996 Urine collection bags

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Connector
3	Hanger
4	Bag

Urine Bags (No Outlet)

STERILE

RE-ORDER	DESCRIPTION	TUBING LENGTH	BAG CAPACITY	CARTON/QTY
UR062004	Urine Bag NRV, with Hanger (no outlet)	120cm	2000mL	50pcs