

Urine Bag No Outlet



PRODUCT DATA SHEET

DESCRIPTION

Our basic Urine Bags are used to collect and contain urine via a catheter or leg bag and made from durable material. They also feature clear volume graduation to aid simple volume measurement. They are available with hanger (as pictured above) or without hanger .

FEATURES AND BENEFITS

- Durable material
- Clear volume graduation
- With non-return valve
- With or without hanger
- No outlet

PRODUCT DATA

- Tube length: 120cm
- Tube ID: 6.5mm
- Tube OD: 9mm
- 2000mL bag
- Material: PVC, PP, PE
- DEHP: Yes
- Sterile: EO
- Shelf life: 5 Years
- ARTG Number (no outlet): 271114

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
UR062004	Urine Bag - NRV Hanger (no outlet)	120 cm Tubing	50
UR062005	Urine Bag - NRV (no outlet)	120cm Tubing	250