

Urine Bag Hanger

PRODUCT DATA SHEET



DESCRIPTION

Our Urine Bag Hanger is designed to simplify the attachment of urine bags to beds. Applying the urine bag can easily be achieved by using the two prongs on the urine hanger to attach to the corresponding perforated holes on the urine bag. The bag tubing can then be secured to the hanger using the pre-formed holder that is located between the bed attachment.


FEATURES AND BENEFITS

- Reinforced double hanger for effective securement
- Easy to use prongs to attach the urine bag
- Suitable for all major bags
- Pre-formed securement clasp to secure the urine bag tubing in place

PRODUCT DATA

- Material: PP
- Latex free
- Sterile: No
- DEHP free
- Shelf life: 5 Years

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	CARTON/QTY
UR063003NS	Urine Bag Hanger	50pcs