

Urine Bag Hanger


DESCRIPTION

Our Urine Bag Hanger is designed to simplify the attachment of urine bags to beds. Applying the urine bag can easily be achieved by using the two prongs on the urine hanger to attach to the corresponding perforated holes on the urine bag. The bag tubing can then be secured to the hanger using the pre-formed holder that is located between the bed attachment.

FEATURES & BENEFITS

- Reinforced double hanger for effective securement
- Easy to use prongs to attach the urine bag
- Suitable for all major bags
- Pre-formed securement clasp to secure the urine bag tubing in place

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Preformed holder

PRODUCT DATA

- Sterile: No
- Not Made with DEHP
- Not Made with Natural Rubber Latex
- Material: PP
- Shelf life: 5 years

ISO STANDARDS

- ISO 13485:2016 - Quality management systems
- EN ISO 15223-1:2016 Symbols for use in the labelling of medical devices
- EN ISO 14644-1 :2015 Cleanrooms and Associated Controlled Environments - Part 1: Classification of Air Cleanliness ISO 14644-1:1999

Urine Bag Hanger

STERILE

RE-ORDER

DESCRIPTION

CARTON/QTY

UR063003NS

Urine Bag Hanger

50