

Open Ended 2-Way Foley Catheters

DESCRIPTION

With its open ended design and bilateral eyelets, this 2-Way Foley Catheter is designed to deliver the high flow rates for urine drainage necessary to support patients with a strong tendency for catheter blockages.

FEATURES & BENEFITS

- Intended for Urethral and SPC use
- Suitable for male and female patients
- Symmetrical balloon
- Colour coded
- Radiopaque marking
- Long indwelling time of eight weeks
- Translucent silicone


PRODUCT DATA

- Length: 40cm
- Side holes: 2
- Balloon size: 10mL
- Not Made with Natural Rubber Latex
- Not Made with DEHP
- Shelf life: 3 years
- ARTG number: 211796
- Sterile EO

ISO STANDARDS

- ENISO 14155:2003 Clinical investigation for medical devices for human subjects
- ISO 10993-10:2002 Biological evaluation of the medical devices-Part 10: Tests for irritation and sensitisation
- EN550:1994 Sterilisation of medical devices-Validation and routine control of ethylene oxide sterilisation

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Open Ended
2	Side Eyelet
3	Balloon
4	X-Ray Line
5	Shaft
6	Funnel
7	Colour Code
8	Check Valve

Open Ended 2-Way Foley Catheters

Open Ended 2-Way Foley Catheter						STERILE
RE-ORDER	DESCRIPTION	TUBING LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
UR074014	Open Ended 2-Way Foley Catheter	40cm	10mL	14Fr	Green	10
UR074016	Open Ended 2-Way Foley Catheter	40cm	10mL	16Fr	Orange	10
UR074018	Open Ended 2-Way Foley Catheter	40cm	10mL	18Fr	Red	10
UR074020	Open Ended 2-Way Foley Catheter	40cm	10mL	20Fr	Yellow	10
UR074022	Open Ended 2-Way Foley Catheter	40cm	10mL	22Fr	Purple	10
UR074024	Open Ended 2-Way Foley Catheter	40cm	10mL	24Fr	Blue	10