


# Open Ended Silicone 2-Way Foley Catheters

## PRODUCT DATA SHEET



### DESCRIPTION

With its open ended design and bilateral eyelets, this 2-Way Foley Catheter is designed to deliver the high flow rates for urine drainage necessary to support patients with a strong tendency for catheter blockages.

### FEATURES AND BENEFITS

- Intended for Urethral and SPC use
- Suitable for male and female patients
- Symmetrical balloon
- Colour coded
- Radiopaque marking
- Long indwelling time of eight weeks
- Translucent silicone

### PRODUCT DATA

- Length: 40cm
- Side holes: 2
- Balloon size: 10mL
- Material: Silicone
- DEHP free
- Latex free
- Shelf life: 3 years
- ARTG Number: 211796

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	GAUGE	BALLOON CAPACITY	COLOUR	CARTON/QTY
UR074014	Open Ended 2-Way Foley Catheter	40cm	14Fr	10mL	Green	10pcs
UR074016	Open Ended 2-Way Foley Catheter	40cm	16Fr	10mL	Orange	10pcs
UR074018	Open Ended 2-Way Foley Catheter	40cm	18Fr	10mL	Red	10pcs
UR074020	Open Ended 2-Way Foley Catheter	40cm	20Fr	10mL	Yellow	10pcs
UR074022	Open Ended 2-Way Foley Catheter	40cm	22Fr	10mL	Purple	10pcs
UR074024	Open Ended 2-Way Foley Catheter	40cm	24Fr	10mL	Blue	10pcs