

# Leg Bag

## T-Tap with Silicone Lined Straps

### DESCRIPTION

The Leg Bags are designed to provide discreet storage of urine under clothing. They are attached to the leg with silicone lined straps to prevent the bag from sliding down the leg. The bags are made from non-woven material, soft to the skin.

Available in; 350mL, 500mL and 750mL versions with 10cm, 30cm or 50cm tube with bonded stepped connector, non-return valve and T-Tap.

### FEATURES & BENEFITS

- Discreet underneath clothing
- Easy to use T-Tap to empty bag
- Available in 350mL, 500mL & 750mL bag sizes
- Bonded stepped connector to tubing
- Non-return valve
- Tubing lengths: 10cm, 30cm & 50cm
- Comfortable silicone lined straps for non-slipping
- Soft non-woven backing


### PRODUCT DATA

- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Material: PVC, PP, PE & Silicone
- Straps: Silicone lined strap with fabric fastener
- Sterile: EO
- Shelf life: 5 years
- ARTG Number: 267876

### ISO STANDARDS

- ISO 8669-2 Urine collection bags – Part 2: Determination of dimensions
- ISO 10993 Biological evaluation of medical device
- ISO 11135 Sterilisation of health care products - Ethylene oxide

### PRODUCT DIAGRAM



ITEM	DESCRIPTION	5	Outlet Tubing
1	Protective Cap	6	T-valve
2	Connector	7	Non Return Film
3	Inlet Tubing	8	Silicone Tube
4	Bag Body	9	Silicone Strap

# Leg Bag

## T-Tap with Silicone Lined Straps

Leg Bags				STERILE
RE-ORDER	DESCRIPTION	LENGTH (CM)	COLOUR	CARTON/QTY
<b>UR135010</b>	LEG BAG - T-Tap NRV Short Tube with Silicone Lined Straps	10cm	350mL	80
<b>UR135030</b>	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	30cm	350mL	80
<b>UR135050</b>	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	50cm	350mL	80
<b>UR150010</b>	LEG BAG - T-Tap NRV Short Tube with Silicone Lined Straps	10cm	500mL	80
<b>UR150030</b>	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	30cm	500mL	80
<b>UR150050</b>	LEG BAG - T-Tap NRV Long Tube with Silicone Lined Straps	50cm	500mL	80
<b>UR175010</b>	LEG BAG - T-Tap NRV Short Tube with Silicone lined Straps	10cm	750mL	80
<b>UR175030</b>	LEG BAG - T-Tap NRV Long Tube with Silicone lined Straps	30cm	750mL	80
<b>UR175050</b>	LEG BAG - T-Tap NRV Long Tube with Silicone lined Straps	50cm	750mL	80